





How to: create a Vocabulary Set 3: adding content to multiple cells (option 1)

Quickly add content to lots of cells by using the symbol importer.




Tap here to add content by selecting symbols.




Search by typing a word or phrase.


Or browse through different topics.




Select the symbol or symbols needed and tap **Done**.



Repeat the procedure until you have added all the symbols required.



To move a cell around, tap and hold on the cell and drag it.



Tap **Select** to choose multiple cells.

Once selected change their colour.

Your Vocabulary Set is now ready to use.

You can select symbols in the symbol selection popover.